

# VENOCATH - SMI

## INTRAVENOUS CANNULA CATHETER



- ☑ Catheter made of resistant FEP (Fluorinated ethylene propylene), thin-walled with optimal flow rates.
- ☑ Lengths specially adapted for veterinary use.
- ☑ This FEP comply with USP Biocompatibility requirements.
- ☑ Recommended catheter dwell time up to 3 days.

SMI REF	SPECIFICATION	G / DIA / L	COLOUR	FLOW ML/MIN	PC./BOX
ZAG26.06.18	G26 0,60 x 19 mm	26 / 06 / 19	violet	17 ml	100
ZAG24.07.19	G24 0,70 x 19 mm	24 / 07 / 19	yellow	23 ml	100
ZAG22.09.25	G22 0,90 x 25 mm	22 / 09 / 25	blue	36 ml	100
ZAG20.11.25	G20 1,10 x 25 mm	20 / 11 / 25	pink	59 ml	100
ZAG20.11.32	G20 1,10 x 32 mm	20 / 11 / 32	pink	56 ml	100
ZAG18.13.32	G18 1,30 x 32 mm	18 / 13 / 32	green	94 ml	100
ZAG18.13.45	G18 1,30 x 45 mm	18 / 13 / 45	green	90 ml	100
ZAG17.15.45	G17 1,50 x 45 mm	17 / 15 / 45	white	142 ml	100
ZAG16.17.45	G16 1,70 x 45 mm	16 / 17 / 45	grey	200 ml	100
ZAG14.20.45	G14 2,00 x 45 mm	14 / 20 / 45	orange	305 ml	100
ZAG14.20.70	G14 2,00 x 70 mm	14 / 20 / 70	orange	290 ml	50
ZAG13.23.50	G13 2,35 x 50 mm	13 / 23 / 50	pink	330 ml	50
ZAG13.23.80	G13 2,35 x 80 mm	14 / 23 / 80	pink	315 ml	50
ZAG12.27.50	G12 2,70 x 50 mm	12 / 27 / 50	blue	420 ml	50
ZAG12.27.80	G12 2,70 x 80 mm	12 / 27 / 80	blue	375 ml	50

- ☑ AISI 304 stainless steel cannula.
- ☑ Catheter connection: Luer lock connector.
- ☑ Double grip plates for secure guidance.
- ☑ Flexible wings to secure fixation.
- ☑ Colour coding for different sizes.
- ☑ Hydrophobic filter for fast reflux. Transparent flash-back chamber.
- ☑ Two radio opaque lines.

For quality  
&  
safety

**smi**

[www.sutures.be](http://www.sutures.be)